




## Small Volumes Made Easy

The VersaPump™ Infuser is a reusable constant pressure mechanical syringe pump controlled by a proprietary variable rate flow device.



## Highlights

The VersaPump™ is a reusable device that offers the clinician and patient a way to infuse while at home, promoting comfort, ease of use, and mobility



## Quick Reference

See instructions for use for complete guide.



**OPEN** the pump to prepare for syringe loading



**LOAD** either a standard or pre filled syringe making sure the syringe flanges are securely locked in place.



**CLOSE** to load the infuser and prepare for infusion. The infuser is validated for a lifetime of use.



**BEGIN** after priming and infusion site preparation per the steps in the IFU. Begin your infusion and leverage the carrying case as needed to enhance mobility.

## Features

Portable, lightweight and quiet

Mechanically powered

Designed for use with the proprietary VersaRate® Plus flow controller or Infuset™ fixed rate sets

## Benefits

Pump easily opens and closes vertically, with no twisting of a knob

Simplicity for clinicians to train patients for home use

Made specifically for lower volume infusions

## Advantages

Offers control of infusion therapy without the cost of an electronic pump

VersaPump™, when used in conjunction with the VersaRate® Plus, gives patients flexibility, and can replace rate control sets

Flow rate charts are available for specific applications

## Ordering Information

For orders and inquires, contact your local EMED representative or EMED's customer service department at 916-932-0071. Orders can also be placed via email at [orders@emedtc.com](mailto:orders@emedtc.com).

EMED Reorder Number

Item Description

Qty Per Box

FP-0010209

VersaPump™ Mechanical Infusion Pump (20/30ml)

1